

Example of ELISA

Example-1 Detection of HBsAgS by sandwich ELISA


Solid phase Ab: anti-HBsAg Ab, Rabbit polyclonal)
 Antigen sample: HBsAg-XT (Cat #BCL-AG-002)
 1st Ab(anti-HBsAg Ab, Mouse monoclonal (Diluted by **Soln. A**)
 or anti-HBsAg Ab, Mouse monoclonal HRP (Diluted by **Soln. B**)
 2nd Ab: anti-mouse IgG, goat polyclonal, HRP (Diluted by **Soln. B**)

Example -2 Detection of anti- HBsAgAb by sandwich ELISA


Solid phase antigen (HBsAg, HBsAg-XT, Cat #BCL-AG-002)
 Ab sample: Anti-HBsAg Ab, mouse mono. (Diluted by **Soln. A**)
 ALP-labeled antigen: AP-BNC-ZZ (Diluted by **Soln. B**)

Results of Example-1

1st and 2nd Ab used



1st Ab alone used



Result of example-2

